

DISPOSABLE FOAM ELECTRODES WITH OFF-CENTERED SNAP CONNECTION 50 mm - round - adult



Product Code	33309
Unit of sale	box of 1000 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

RDM (NSIS)	2262590
CND	C020501
EAN/UPC	8054953341224
GMDN	35035

Description

Disposable electrodes with off-centered snap connector

FIAB code: F2017
Material: PE-Foam
Size: 50 mm
Gel: Liquid
Patient: Adult
Minimun order: box of 1,000 (20 bags of 50)

Application: short and medium term monitoring, stress ECG recording and holter

All electrodes have silver chloride sensor