

**DISPOSABLE PE-FOAM ELECTRODES 28-44 mm -
rectangular - adult/ped.**



Product Code	33311
Unit of sale	box of 2000 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

RDM (NSIS)	2254828
CND	C020501
EAN/UPC	8054953342191
GMDN	35035

Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F9046
- Material: PE-Foam
- Size: 28x44 mm
- Gel: Solid
- Shape: Rectang.
- Connection: Snap
- Patient: Adult, Pediatric
- Application: rest, stress, ECG holter
- Minimum order: Box of 2,000 pcs (20 bags of 100)