

DISPOSABLE RADIOTRANSSPARENT PE-FOAM ELECTRODES 28-44 mm - rectangular - adult/ped.



| | |
|---------------|-----------------|
| Product Code | 33312 |
| Unit of sale | box of 2000 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
|--------|-----|

| | |
|------------|---------------|
| RDM (NSIS) | 2262585 |
| CND | C020501 |
| EAN/UPC | 8054953343495 |
| GMDN | 35035 |

Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F9046 RM
- Radiolucent and amagnetic
- Material: PE-Foam
- Size: 28x44 mm
- Gel: Solid
- Shape: Rectang.
- Connection: Carbon Snap
- Patient: Adult, Pediatric
- Application: rest, stress, ECG holter, MRI
- Minimum order: Box of 2,000 pcs (20 bags of 100)