



Product Code	33313
Unit of sale	box of 2000 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2251899
CND	C020501
EAN/UPC	8054242449655
PARAF	962639047
GMDN	35035

Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F9089P-100
- Material: Foam
- Size: Ø 32 x 36 mm
- Gel: Solid
- Shape: Oval
- Connection: Snap
- Patient: Pediatric
- Application: rest, stress, ECG holter
- Minimum order: Box of 2,000 pcs (20 bags of 100)