



Product Code	33314
Unit of sale	box of 2000 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2251898
CND	C020501
EAN/UPC	8054242449648
PARAF	987474956
GMDN	35035

Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F9089-100
- Material: Foam
- Size: Ø 36 x 40 mm
- Gel: Solid
- Shape: Oval
- Connection: Snap
- Patient: Adult
- Application: rest, stress, ECG holter
- Minimum order: Box of 2,000 pcs (20 bags of 100)