

DISPOSABLE FOAM ELECTRODES 42-45 mm - adult - liquid gel



Product Code	33315
--------------	-------

Unit of sale	box of 1250 pcs
--------------	-----------------

Minimum order	1
---------------	---

Type	Medical device
------	----------------

Class	I
-------	---

UK-REP	No
--------	----

CH-REP	No
--------	----

RDM (NSIS)	2251900
------------	---------

CND	C020501
-----	---------

EAN/UPC	8054953347141
---------	---------------

Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F9089GL
- Material: Foam
- Size: Ø 42x45 mm
- Gel: Liquid
- Shape: Oval
- Connection: Snap
- Patient: Adult
- Application: rest, stress, ECG holter
- Minimum order: box of 1,250 (25 bags of 50)