

# ULTRA CHEST ELECTRODES diameter 30 mm



Product Code	33324
Unit of sale	set of 6 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2348589
CND	C020501
EAN/UPC	8033003261569
UMDNS	11-439
GMDN	42489

## Description

### UNIVERSAL REUSABLE ECG ELECTRODES

Set of 6 chest electr. Ultra Ø 30 mm - silver coated

PVC blue balls Ø 33 mm  
Electrode Ø 4 mm  
4 mm elastic connector