



Product Code	33327
Unit of sale	box of 1200 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2251933
CND	C020501
EAN/UPC	8033695725455
PARAF	987474968
GMDN	35035

Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F2080
- Material: Non woven
- Size: Ø 55 mm
- Gel: Solid
- Shape: Round
- Connection: Snap
- Patient: Adult
- Application: rest, long term, stress, ECG holter
- Minimum order: Box of 1,200 pcs (24 bags of 50)