



Product Code	33341
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2207687
CND	C020501
PARAF	977813854
GTIN	38033695726514
GMDN	35035

## Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F3001ECG
- Material: Gelled
- Size: 34 x 23 mm
- Gel: Solid
- Shape: TAB
- Connection: Alligator
- Patient: Adult
- Application: rest ECG
- Minimum order: Bag of 100 pcs.