



| | |
|---------------|----------------|
| Product Code | 33343 |
| Unit of sale | box of 520 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| RDM (NSIS) | 2251910 |
| CND | C020501 |
| EAN/UPC | 8033695726148 |
| GMDN | 35035 |

Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F3001NECG
- Material: Gelled
- Size: 34 x 13 mm
- Gel: Solid
- Shape: TAB
- Connection: Alligator
- Patient: Pediatric, neonatal
- Application: rest ECG
- Minimum order: Box of 520 pcs (13 bags of 40)