



Product Code	33344
Unit of sale	box of 1200 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2262588
CND	C020501
EAN/UPC	8033695726872
PARAF	987474943
GMDN	35035

Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F9067
- Material: PE - Foam
- Size: Ø 48 x 50 mm
- Gel: Liquid
- Shape: Oval
- Connection: Snap
- Patient : adulte
- Application : repos, stress, holter ECG
- Minimum order: Box of 1,200 pcs (24 bags of 50)