

TAB ELECTRODES - pre-gelled, disposable (50 bags of 100)



Product Code	33346
Unit of sale	box of 5000 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2207687
CND	C020501
EAN/UPC	8033695726513
GMDN	35035

Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F3001ECG
- Material: Gelled
- Size: 34 x 23 mm
- Gel: Solid
- Shape: TAB
- Connection: Alligator
- Patient: Adult
- Application: rest ECG
- Minimum order: Box of 5,000 pcs (50 bags of 100)