

# LIMB CLAMP ELECTRODES paediatric



Product Code	33363
Unit of sale	set of 4 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2243337
CND	C020501
UMDNS	11-439
PARAF	987471861
GTIN	08033003263525
GMDN	42489

## Description

### UNIVERSAL REUSABLE ECG ELECTRODES

Set of 4 limb clamp electrodes (4 colours) paediatric.

20 x 90 mm paediatric clamp

Universal connection: hole with screw (4mm) and clip.

Silver coated plastic material conform to EN 60601-2-25