

LIMB CLAMP ELECTRODES paediatric



Product Code 33363

Unit of sale set of 4 pcs

Minimum order 1

Type Medical device

Class I

RDM (NSIS) 2243337

CND C020501

UMDNS 11-439

PARAF 987471861

GMDN 42489

Description

UNIVERSAL REUSABLE ECG ELECTRODES

Set of 4 limb clamp electrodes (4 colours) paediatric.

20 x 90 mm paediatric clamp

Universal connection: hole with screw (4mm) and clip.

Silver coated plastic material conform to EN 60601-2-25