

CHEST ELECTRODES diameter 24 mm



Product Code	33365
Unit of sale	set of 6 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2501259
CND	C020501
EAN/UPC	8023279333657
UMDNS	11-439
PARAF	903177970
GTIN	08023279333657
GMDN	42489

Description

UNIVERSAL REUSABLE ECG ELECTRODES

Set of 6 chest electrodes Ø 24 mm - silver coated

PVC blue balls Ø 33 mm

4 mm metal connector hole with screw and clips