

DISPOSABLE PE-FOAM ELECTRODES 40 mm - round - pediatric



Product Code	33370
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2251912
CND	C020501
PARAF	974039962
GTIN	08033695726094
GMDN	35035

Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F9060P
- Material: PE - Foam
- Size: Ø 40 mm
- Gel: Solid
- Shape: Round
- Connection: Snap
- Patient: Pediatric
- Application: rest, long term, stress, ECG holter
- Minimum order: Bag of 50 pcs.