

# DISPOSABLE PE-FOAM ELECTRODES 48-50 mm - oval - adult



Product Code	33371
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2251911
CND	C020501
PARAF	974039950
GTIN	08033695724304
GMDN	35035

## Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F9060
- Material: PE - Foam
- Size: Ø 48 x 50 mm
- Gel: Solid
- Shape: Oval
- Connection: Snap
- Patient: Adult
- Application: rest, long term, stress, ECG holter
- Minimum order: Bag of 50 pcs.