

# DISPOSABLE PE-FOAM ELECTRODES 28-44 mm - rectangular - adult



Product Code	33372
Unit of sale	box of 1000 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2251903
CND	C020501
EAN/UPC	8054242440096
GMDN	35035

## Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F9047/4FL
- Material: PE - Foam
- Size: 28 x 44 mm
- Gel: Liquid
- Shape: Rectangular
- Connection: Plug 4 mm
- Patient: Adult
- Application: rest ECG
- Minimum order: Box of 1,000 pcs (20 bags of 50)