

DISPOSABLE FOAM ELECTRODES - oval 23-30 mm with wire 35 cm - pediatric



| | |
|---------------|----------------|
| Product Code | 33373 |
| Unit of sale | box of 150 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|---------------|
| RDM (NSIS) | 2254711 |
| CND | C020501 |
| EAN/UPC | 8033003261682 |
| GMDN | 35035 |

Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F9053N
- Material: Foam
- Size: Ø 23 x 30 mm
- Gel: Solid
- Shape: Oval
- Connection: 4 mm socket
- Patient: Pediatric, neonatal
- Application: monitoring ECG recording
- Minimum order: Box of 1,50 pcs (5 bags of 30)