

# DISPOSABLE FOAM ELECTRODES - oval 23-30 mm with wire 50 cm - pediatric



Product Code	33377
Unit of sale	box of 150 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2254769
CND	C020501
EAN/UPC	8033003260357
GMDN	35035

## Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F9058N
- Material: Foam
- Size: Ø 23 x 30 mm
- Gel: Solid
- Shape: Oval
- Connection: 1.5 mm socket
- Patient: Pediatric, neonatal
- Application: long monitoring ECG recording
- Minimum order: Box of 1,50 pcs (5 bags of 30)