



Product Code	33343
Unit of sale	box of 520 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2251910
CND	C020501
EAN/UPC	8033695726148
GMDN	35035

Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F3001NECG
- Material: Gelled
- Size: 34 x 13 mm
- Gel: Solid
- Shape: TAB
- Connection: Alligator
- Patient: Pediatric, neonatal
- Application: rest ECG
- Minimum order: Box of 520 pcs (13 bags of 40)