

NOSE CLIP - rubber - large



Product Code	33477
Unit of sale	box of 100 pcs
Minimum order	1
Type	No medical device
Class	

EAN/UPC	8023279334777
---------	---------------

Description

NOSE CLIPS - rubber large - 70x45 mm - bag of 100

Gima rubber or foam-padded disposable nose clips help prevent patients from inhaling or exhaling through the nose during spirometry. The tension in the nose clip body keeps it firmly in place while the rubber or foam pads are designed for the best patient comfort. Made of thermoplastic, latex free.

Rubber nose clip is available in 2 sizes: standard or large.

Made in Italy.