



Product Code	33559
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	360711
CND	Z12030599
EAN/UPC	8023279335590
GMDN	17882

Description

Semi automatic external defibrillators designed for public use by minimally trained people. They provide easy and simple indications for rescue operations. Designed to treat Ventricular Fibrillation (VF) and Fast Ventricular Tachycardia (VT). Supplied with lithium battery, set of adult pads and manual. According to AHA/ERC Guidelines 2010. Made in Italy.

Other languages available in 20 days: DE, PL, RU (manual and voice), PT, CZ, FI, GR, LT, NL, RO, TR, KR, TH (voice with English manual).

Technical Specifications

- Operating mode: Semi-automatic
- Waveform: BTE (Biphasic Exponential Truncated) impedance compensated
- Energy: 200 J fixed
- Charge control: Automatic by software (Arrhythmia Detection System and charging control)
- Charge time: < 8 seconds, typical with new battery
- Shock delivery vector: adult pads in the anterior (Lead II) position or reduced energy paediatric pads in anterior-posterior position
- Voice instructions: Vocal and optical message guide
- Automatic manual tests: Battery insertion test and daily self tests
- Long life lithium battery: 12 Volt DC, 4.2 Ah, LiMnO₂ disposable
- Battery capacity: 200 shocks or 4 hours of continuous operating time
- Data recording, stored: first 60 minutes of ECG and the entire incident's event
- Size-weight: 30 x 28 x h 9 cm - 2.3 kg with battery