


# COMPATIBLE PADS for defibrillator Agilent-Philips see also 55020



Product Code	33576
Unit of sale	kit of 2 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2468982
CND	C020401
EAN/UPC	8033003260319
PARAF	951450384
GMDN	47055

## Description

### ADULT Compatible pads - disposable

A range of compatible pads BSI certified for patients >8 years (25 kg) for most defibrillators on the market.

Excellent adherence and adequate gel thickness.

Includes a couple of pads, cable and connector.

Maximum shelf life 30 months (24 months for 33584) from production date.

Latex free.

Compliance to regulations: Directive 93/42/EEC and amendments. Class IIb medical device.

ANSI/AAMI: DF-80; IEC/EN: 60601-1, 60601-2-4; ISO: 10993-1.

Multilingual envelope: GB, FR, IT, ES, PT, DE, NL, SE, GR, RU.

Made in Italy.

### Compatible with:

- Agilent-Philips
    - Heartstream/HeartStart XL M4735A
    - Heartstream/HeartStart XLT M3500B
    - HeartStart XL 861290
    - HeartStart 4000 M5500B
- (all with Philips cable M3508A)