

COMPATIBLE PADS for defibrillator Agilent, H-P, Laerdal see also 55022



Product Code	33577
Unit of sale	kit of 2 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2468987
CND	C020401
EAN/UPC	8033003260340
GMDN	47055

Description

ADULT Compatible pads - disposable

A range of compatible pads BSI certified for patients >8 years (25 kg) for most defibrillators on the market. Excellent adherence and adequate gel thickness.

Includes a couple of pads, cable and connector.

Maximum shelf life 30 months (24 months for 33584) from production date.

Latex free.

Compliance to regulations: Directive 93/42/EEC and amendments. Class IIb medical device.

ANSI/AAMI: DF-80; IEC/EN: 60601-1, 60601-2-4; ISO: 10993-1.

Multilingual envelope: GB, FR, IT, ES, PT, DE, NL, SE, GR, RU.

Made in Italy.

Compatible with:

- Agilent Hewlett-Packard
 - Codemaster 100-XE-XL-XL+
 - Heartstream XL, XLT, XL+

- Laerdal
 - HeartStart 4000