


COMPATIBLE PADS for defibrillator Cardiac Science, GE see also 55026



Product Code	33580
Unit of sale	kit of 2 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2468894
CND	C020401
EAN/UPC	8033695725196
PARAF	988194686
GMDN	47005

Description

ADULT Compatible pads - disposable

A range of compatible pads BSI certified for patients >8 years (25 kg) for most defibrillators on the market.

Excellent adherence and adequate gel thickness.

Includes a couple of pads, cable and connector.

Maximum shelf life 30 months (24 months for 33584) from production date.

Latex free.

Compliance to regulations: Directive 93/42/EEC and amendments. Class IIb medical device.

ANSI/AAMI: DF-80; IEC/EN: 60601-1, 60601-2-4; ISO: 10993-1.

Multilingual envelope: GB, FR, IT, ES, PT, DE, NL, SE, GR, RU.

Made in Italy.

Compatible with:

- Cardiac Science
 - Powerheart AED 9200 series
 - Powerheart AED G3 9300 series (G3, G3 Plus, G3 Pro)

- GE
 - Responder AED - AED Pro