


COMPATIBLE PADS for defibrillator Mediana, Tecno-Gaz see also 55046



Product Code	33584
Unit of sale	kit of 2 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2474133
CND	C020401
EAN/UPC	8054242449587
GMDN	47005

Description

ADULT/Pediatric Compatible pads - disposable
(CPR guidelines state 8 years as limit age for pediatric treatment, corresponding to about 25 kg according to standard growth curves)

A range of compatible pads BSI certified for patients >8 years (25 kg) for most defibrillators on the market.

Excellent adherence and adequate gel thickness.

Includes a couple of pads, cable and connector.

Maximum shelf life 30 months (24 months for 33584) from production date.

Latex free.

Compliance to regulations: Directive 93/42/EEC and amendments. Class IIb medical device.

ANSI/AAMI: DF-80; IEC/EN: 60601-1, 60601-2-4; ISO: 10993-1.

Multilingual envelope: GB, FR, IT, ES, PT, DE, NL, SE, GR, RU.

Made in Italy.

Compatible with:

- Mediana
- HeartOn AED A15 (patient >10 kg)
- Tecno-Gaz
- Tecnoheart Plus (patient >10 kg)