


# COMPATIBLE PADS for defibrillator Medtronic, Osatu Bexen see also 55048



|               |                |
|---------------|----------------|
| Product Code  | 33587          |
| Unit of sale  | kit of 2 pcs   |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II B           |

|            |               |
|------------|---------------|
| UK-REP     | Yes           |
| CH-REP     | Yes           |
| RDM (NSIS) | 2468985       |
| CND        | C020401       |
| EAN/UPC    | 8033003260333 |
| GMDN       | 47055         |

## Description

ADULT Compatible pads - disposable

A range of compatible pads BSI certified for patients >8 years (25 kg) for most defibrillators on the market.  
Excellent adherence and adequate gel thickness.  
Includes a couple of pads, cable and connector.  
Maximum shelf life 30 months (24 months for 33584) from production date.  
Latex free.  
Compliance to regulations: Directive 93/42/EEC and amendments. Class IIb medical device.  
ANSI/AAMI: DF-80; IEC/EN: 60601-1, 60601-2-4; ISO: 10993-1.  
Multilingual envelope: GB, FR, IT, ES, PT, DE, NL, SE, GR, RU.  
Made in Italy.

Compatible with:

- Bexen Cardio  
- Reanibex 200 - 300 - 500 - 700 - 800
- Cardioline  
- Elife 700
- Mindray  
- Beneheart D1 - D3 - D6
- Physiocontrol  
- Lifepak 9 connector Quick-Combo  
- Lifepak 10 - 12 - 15 - 20 - 20E - 500 AED pink or grey socket - 1000 AED