



Product Code	33613
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2351778
CND	Z121401
EAN/UPC	8054187380341

## Description

Piccolo Basic screening audiometer - air + mask PC for Pc use

Piccolo is a basic diagnostic audiometer offering air and bone audiometry functions. Piccolo can be controlled from a computer running Windows. Lightweight, compact design, powerful, userfriendly software and Noah connectivity combine to make Piccolo the ideal choice for professionals on the move.

- computer control through the Maestro software.

Compatible for Windows 7, 8, 10

- small and light weight

- power supply through the USB connection, when controlled by computer.

Software language: GB, FR, IT, ES, DE, PL, PT, RU.

CD manual and quick guide GB, FR, IT, ES, DE, PL, PT, NL.

Piccolo needs annual calibration by authorized centre. Please check on [www.gimaitaly.com](http://www.gimaitaly.com) list of countries where you can find a centre.

Made in Italy.

## Technical Specifications

Type of device: Pc control

Class EN 60645-1, ANSI S3.6: 4-screening

Output: AC

Masking: NBN

Warble tone: Yes

Pulsed tone: 0.5 Hz, 1 Hz, 2 Hz

Frequency range: 125 - 8,000 Hz

Maximum output intensities: AC: -10 to 120 dB HL

Attenuator step: 5 dB

Automatic audiometry: HW auto threshold

Talk over: External microphone

Connection to PC: USB [driverless]

Controlling software: Maestro module

Dimensions: 16x16x3cm (16.3x6.3x1.2 in)

Weight: 300 g (10.6 oz)

Operating voltage: 100-240 V, 50/60 Hz

## Standard accessories

- DD445 supra-aural headphones
- Patient response switch

- Clip-on microphone for operator to-patient communication
- Inventis Daisy Patient Record Manager, with Daisy Maestro module
- Noah Maestro module
- USB connection cable
- Medical grande power supply
- User Manual