

PICCOLO SPEECH DIAGNOSTIC AUDIOMETER - air+bone+mask



Product Code	33618
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2351780
CND	Z121401
EAN/UPC	8054187380389
GMDN	37503

Description

Piccolo Speech diagnostic audiometer - air + bone + mask

It includes Master Hearing Aid

Piccolo is a basic diagnostic audiometer offering air, bone and speech audiometry functions.

Piccolo can be controlled from a computer running Windows.

Lightweight, compact design, powerful, user-friendly software and Noah connectivity combine to make Piccolo the ideal choice for professionals on the move.

Computer control through the Maestro software.

Compatible for Windows 7, 8, 10

Power supply through the USB connection, when controlled by computer.

Software language: GB, FR, IT, ES, DE, PL, PT, RU.

CD manual and quick guide GB, FR, IT, ES, DE, PL, PT, NL.

Piccolo needs annual calibration at any audiometric calibration centre.

Made in Italy.

Technical Specifications

Type of device: PC control

Class IEC 60645-1/ANSI S3.6: 3B diagnostic

Output: AC/BC, Live Speech

Input: - - EXT 1, EXT 2, MIC

Masking: NBN, WN, SN

Warble tone: Yes

Pulsed tone: 0.5 Hz, 1 Hz, 2 Hz

Frequency range: 125 - 8,000 Hz

Maximum output intensities: AC: -10 to 120 dB HL / BC: -10 to 80 dB HL

Attenuator step: 1 dB, 3 dB, 5 dB

Automatic audiometry: HW auto threshold

Talk over: External microphone

Connection to PC: USB (driverless)

Compatible software: Maestro module

Dimensions: 16x16x3cm (6.3x6.3x1.2 in)

Weight: 300 g (10.6 oz)

Operating voltage: 100-240 V, 50/60 Hz