

TRIANGLE BASIC SCREENING AUDIOMETER - air



Product Code	33619
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2351781
CND	Z121401
EAN/UPC	8054187380778

Description

TRIANGLE BASIC SCREENING AUDIOMETER - air

Triangle is a battery-operated, portable screening audiometer, capable of performing automatic and manual pure tone audiometry exams in air conduction.

Designed to deliver fast, accurate, and reliable pure tone audiometry for healthcare professionals on the move.

Its unique combination of advanced technology and userfriendly design, featuring a capacitive colour touchscreen, allow for intuitive operation, making hearing screening simpler and more efficient than ever.

- Pure Tone Audiometry in AC

- compact and lightweight

- rechargeable lithium battery delivers over 12 hours of continuous use

- internal memory for 100 patients

Made in Italy.

Technical Specifications

Type of device: screening audiometer

Class IEC 60645-1/ANSI S3.6: type 4

Output: AC

Warble tone: Yes

Pulsed tone: 1 Hz, 2 Hz

Frequency range: 125-8,000 Hz

Maximum output intensities: AC: -10 to 100 dB HL

Attenuator step: 5 dB

Automatic audiometry: available algorithms: Hughson Westlake, Quick Search, Fixed Intensity

Talk over: built-in microphone

Compatible software: Maestro

Memory: up to 100 patients

Display: TFT LCD touchscreen display, 240x320 resolution

Display size: 2.8" , 43.2x57.6 mm

Power supply: built-in lithium rechargeable battery

Operating voltage: 100-240 V, 50/60 Hz

Size -weight: 21x16x4 cm - 425 g

Standard accessories

DD45 supra-aural headphones

USB connection cable

Patient response switch

Soft carrying case

Audiogram cards

Medical grade power supply

User manual: GB, FR, IT, ES, DE.

On-line manuals PT, NL, DK, SE, NO, LT, LV, EE, PL, HU, RO, BG, GR, TR available on gimaitaly.com