



Product Code	33643
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2224551
CND	Z12149002
EAN/UPC	8023279336436

Description

PRO 28 BASIC - 96x96x197 cm

Soundproof booth with innovative internal isolating design and acoustic technology for a higher soundproofing performance in the most compact solution.

The internal finishing is composed by a special soundproofing fiber covered with a sanitizable, fireproof, highly resistant acoustic material which provides the booth with a perfect acoustics.

Recommended for hearing screening and diagnostic tests for those who have limited space available.

Made in Italy.

Technical Specifications

External size (cm): 96x96x 197

Internal size (cm): 81x81xh 180

Door size cm: 55xh 168

opening: Left/Right

Windows size cm: 52xh 52

material: 2 - crystal insulating safety glass

Wall thickness cm: 7.5

Weight (kg): 233

Ventilation: no

Tested noise reduction dB at 1,000 Hz: 36.2

Internal illumination power V-Hz: 230 - 50/60

Norms: Class I Medical Devices European Directive EU 2017/745

Standard accessories

- module for silent audiometer cable connection
- internal illumination
- ISPM-15 wooden box for transport
- manual: GB, IT