

PRO 28 INTERMEDIATE AUDIOMETRIC BOOTH 96x96x197 cm with ventilation



| | |
|---------------|----------------|
| Product Code | 33644 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| RDM (NSIS) | 2224551 |
| CND | Z12149002 |
| EAN/UPC | 8023279336443 |

Description

PRO 28 INTERMEDIATE - 96x96x197 cm

Soundproof booth with innovative internal isolating design and acoustic technology for a higher soundproofing performance in the most compact solution.

The internal finishing is composed by a special soundproofing fiber covered with a sanitizable, fireproof, highly resistant acoustic material which provides the booth with a perfect acoustics.

Integrated silent forced ventilation system

Recommended for hearing screening and diagnostic tests for those who have limited space available.

Made in Italy.

Technical Specifications

External size (cm): 96x96x197
Internal size (cm): 81x81x173.5
Door size cm: 55x168
opening: Left/Right
Windows size cm: 52x52
material: 2 - crystal insulating safety glass
Wall thickness cm: 7.5
Weight (kg): 233
Ventilation: yes
Tested noise reduction dB at 1,000 Hz: 36.2
Internal illumination power V-Hz: 230 - 50/60
Norms: Class I Medical Devices European Directive EU 2017/745

Standard accessories

- module for silent audiometer cable connection
- internal illumination
- ventilation system
- ISPM-15 wooden box for transport
- manual: GB, IT