

# PRO 28 PREMIUM AUDIOMETRIC BOOTH 96x96x197 cm with ventilation and window door



- 33651
- tavolino (opzionale)
  - table (optional)
  - table (optionnel)

Product Code	33649
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2224551
CND	Z12149002
EAN/UPC	8023279336498
GMDN	37020

## Description

PRO 28F - 96x96x197 cm

New compact booth with an outstanding level of soundproofing and excellent acoustic thanks to innovative materials.

Internal finishing made of fire proof, washable, sanitizable and very resistant fabric.

Internal ventilation system.

External finishing made of highly resistant, fireproof, washable panel.

Module for cable connection included.

Made in Italy.

Reversible door (left/right) with magnetic closing.

## Technical Specifications

External size (cm): 96x96xh 197

Internal size (cm): 81x81xh 173.5

Door size cm: 55xh 168

opening: Left/Right

Windows size cm: 52xh 52

material: 2 - crystal insulating safety glass

Wall thickness cm: 7.5

Weight (kg): 233

Ventilation: yes

Tested noise reduction dB at 1,000 Hz: 36.2

Internal illumination power V-Hz: 230 - 50/60

Norms: Class I Medical Devices European Directive EU 2017/745

## Standard accessories

- module for silent audiometer cable connection
- internal illumination
- ventilation system
- ISPM-15 wooden box for transport
- manual: GB, IT