



33651
 • tavolino (opzionale)
 • table (optional)
 • table (optionnel)

Product Code	33654
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2227100
CND	Z12149002
EAN/UPC	8023279336542
GMDN	37020

Description

PRO 30 - 107x107x244 cm - Height includes air ventilation system (27 cm)

Interior finishing made of auto extinguishing acoustic pyramid foam in anthracite grey.
 Range with different sizes, with high noise reduction and ventilation system.
 External finishing made of highly resistant, fireproof, washable panel and aluminium frame.
 Module for cable connection included.
 Sanitizing air-exchanging system on request.
 Made in Italy.

Technical Specifications

External size (cm): 107x107xh 244 - Height includes air ventilation system (27 cm)
 Internal size (cm): 92x92xh 201
 Door size cm: 75xh 200
 opening: Right (Left opening available on request)
 Windows size cm: 56xh 88
 material: 2 - crystal insulating safety glass
 Wall thickness cm: 7.5
 Weight (kg): 285
 Ventilation: yes
 Tested noise reduction dB at 1,000 Hz: 34.5
 Internal illumination power V-Hz: 86-265 - 50/60
 Norms: Class I Medical Devices European Directive EU 2017/745

Standard accessories

- module for silent audiometer cable connection
- internal illumination
- ventilation system
- ISPM-15 wooden box for transport
- manual: GB, IT