

SPACER FOR AEROSOL 175 ml - adult mask (5 years +)



Product Code	33672
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2732496
CND	R03010380
EAN/UPC	8023279336726
PARAF	987472685
GMDN	38018

Description

SPACER FOR AEROSOL 175 ml - adult (5 years +)

Large mask recommended for patients 5 years +
Weight: 86 g

Anti-static simple and portable aerosol chamber suitable for people with respiratory diseases.
It can manage asthma, allowing patients to inhale medicines in a more effectively way.
Using a spacer enables to breathe in the mist directly into the lungs and not wasting in the back of the mouth.
It helps providing relief from asthma symptoms.
It shouldn't be used by people with severe asthma or breathing difficulties.
Flow signal whistle to alert patients to breathe slowly included.
Chamber size: 52x131 mm.
Single patient use only, not sterile.
Multilingual box and instructions: GB, FR, IT, ES, PT, DE, NL, DK, SE, LT, PL, CZ, HU, HR, RO, GR.