

# SPACER FOR AEROSOL 175 ml - baby mask (0-18 months)



Product Code	33674
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2732496
CND	R03010380
EAN/UPC	8023279336740
PARAF	987472697
GMDN	38018

## Description

SPACER FOR AEROSOL 175 ml - baby (0-18 months)

Small mask anatomically shaped for a secure seal  
Weight: 63 g

Anti-static simple and portable aerosol chamber suitable for people with respiratory diseases.  
It can manage asthma, allowing patients to inhale medicines in a more effectively way.  
Using a spacer enables to breathe in the mist directly into the lungs and not wasting in the back of the mouth.  
It helps providing relief from asthma symptoms.  
It shouldn't be used by people with severe asthma or breathing difficulties.  
Flow signal whistle to alert patients to breathe slowly included.  
Chamber size: 52x131 mm.  
Single patient use only, not sterile.  
Multilingual box and instructions: GB, FR, IT, ES, PT, DE, NL, DK, SE, LT, PL, CZ, HU, HR, RO, GR.