

CHEST ELECTRODES diam. 15 mm



Product Code	33368
Unit of sale	set of 6 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2501261
CND	C020501
EAN/UPC	8023279333688
UMDNS	11-439
GTIN	08023279333688
GMDN	42489

Description

UNIVERSAL REUSABLE ECG ELECTRODES

Set of 6 chest electrodes Ø 15 mm - silver coated

PVC balls Ø 25 mm

4 mm metal connector hole with screw and clips