

3-LEAD PATIENT CABLE (cable + snap connectors 33748)



Product Code	33711
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

EAN/UPC	8023279337112
---------	---------------

Description

ECG CABLES FOR BM3 (33718, 33719, 33720, 33780, 33781)

3-lead patient cable kit (cable+snap connectors 33748)

- kit composition:
 - Snap connectors (33748)
 - 3-leads cable (33762)
- New model (after 6/2006)