

3-LEAD VET CABLE (cable + 3 clips connectors)



Product Code 33713

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279337136

Description

ECG CABLES FOR BM3 (33718, 33719, 33720, 33780, 33781)

3-lead veterinary cable kit (cable+3 clips connectors)

- kit composition:
 - Vet clip connectors (33760)
 - 3-leads cable (33762)
- New model (after 6/2006)