

METAL CONNECTOR (necessary for mini cuff)



Product Code	33739
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8809276946829
---------	---------------

Description

BM1, BM3, BM5 and BM7 monitors accessories:
Metal connector (necessary for mini cuff codes 32810/11).