


# INFANT PROBE under 3 kg - over 40 kg Nellcor comp.- disposable



Product Code	33755
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	1720266
CND	C900301
EAN/UPC	6945040104346
UMDNS	37808

## Description

Disposable SpO<sub>2</sub> sensor

- For neonate <3 kg (neonate foot)
- For adults > 40 kg (adult finger)

Compatibility:

- NELLCOR (not compatible with Nellcor Oximax)

All sensors are supplied in sterile pack with 70 cm cable.

## Technical Specifications

Range of SpO<sub>2</sub> measurement: 70%~100%  
Accuracy: 70%~100% ±2%; below 70%: unspecified  
Range of pulse measurement: 30~250 bpm  
Accuracy: ±2 bpm of 2% (select larger)  
Optical Sensor:  
red light (wavelength is 650~670 nm, 6.65 mW)  
infrared (wavelength is 880~910 nm, 6.75 mW)