



Product Code	33757
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	1949941
CND	C900301
EAN/UPC	6945040104322
UMDNS	37808
PARAF	987472750

## Description

Disposable SpO2 sensor

- Pediatric 10-50 kg

Compatibility:

- NELLCOR (not compatible with Nellcor Oximax)

All sensors are supplied in sterile pack with 70 cm cable.

## Technical Specifications

Range of SpO2 measurement: 70%~100%  
 Accuracy: 70%~100%  $\pm 2\%$ ; below 70%: unspecified  
 Range of pulse measurement: 30~250 bpm  
 Accuracy:  $\pm 2$  bpm of 2% (select larger)  
 Optical Sensor:  
 red light (wavelength is 650~670 nm, 6.65 mW)  
 infrared (wavelength is 880~910 nm, 6.75 mW)