

DISPOSABLE FOAM ELECTRODES - oval 23-30 mm with  
wire 50 cm - pediatric



|               |                |
|---------------|----------------|
| Product Code  | 33377          |
| Unit of sale  | box of 150 pcs |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2254769       |
| CND        | C020501       |
| EAN/UPC    | 8033003260357 |
| GMDN       | 35035         |

Description

UNIVERSAL DISPOSABLE ECG ELECTRODES - latex free

- FIAB code: F9058N
- Material: Foam
- Size: Ø 23 x 30 mm
- Gel: Solid
- Shape: Oval
- Connection: 1.5 mm socket
- Patient: Pediatric, neonatal
- Application: long monitoring ECG recording
- Minimum order: Box of 1,50 pcs (5 bags of 30)