



Product Code	33796
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	1801762
CND	Z1203020280
EAN/UPC	8809276943439
GMDN	36554

Description

LOFLO SIDESTREAM CO2 SENSOR for intubated and non intubated patients

Flexible and compact CO2 sensor that provides consistent and reliable CO2 monitoring of adult, pediatric and neonatal patients. No calibration required. 50 ml/minute sampling rate.

COMPATIBILITY:

Bionet BM5, BM5 Pro
Creative K15
Mindray Various
Philips Various