

# RESPIRONICS CAPNOSTAT ETCO2 MAINSTREAM SENSOR



Product Code	33797
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2503056
CND	Z1203020480
EAN/UPC	0884838013667
GMDN	36552

## Description

CAPNOSTAT MAINSTREAM CO2 SENSOR for intubated patients

Small, durable and lightweight mainstream CO2 sensor provides accurate and reliable monitoring for all intubated patients from neonates to adults. No calibration required. Fast airway measurement.

### COMPATIBILITY:

Bionet BM5, BM5 Pro  
Creative K15  
Mindray Various  
Philips Various