



Product Code	33871
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2975214
CND	Z11040182
EAN/UPC	6945121407403
UMDNS	15976

## Description

Chison linear probe, center frequency 9.0 MHz for Portable Ultrasound Machines ECO series.

Probe for Chison portable ultrasounds:

- Chison code: L7S-A
- Scanning method: Linear (with 30 mm footprint for superficial tissue exams with higher resolution)
- Main Frequency: 9.0 MHz
- Frequency Range(MHz): 5.3 / 11
- Applications: Small Parts, MSK, Vascular, Pediatric, Orthopaedic