

CHISON 7.5 MHz TRANSRECTAL PROBE for ECO1/2/3/5/6



Product Code	33873
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2975214
CND	Z11040182
EAN/UPC	6945121407410
UMDNS	15976

Description

Probe for Chison portable ultrasounds (For ECO3EXP., ECO5 and ECO6 only)

- Chison code: R7-A
- Scanning method: Trans-rectal
- Main Frequency: 7.5 MHz
- Frequency Range (MHz): 5.0 / 10.0
- Applications: Urology