


CHISON 3.0 MHz MICRO-CONVEX PROBE for ECO1/2/3/5/6



Product Code	33874
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2975214
CND	Z11040182
EAN/UPC	6945121407373
UMDNS	15976

Description

Probe for Chison portable ultrasounds:

- Chison code: MC3-A
- Probe type: Micro-Convex
- Main Frequency: 3.0 MHz
- Frequency Range(MHz): 2.5 / 5.0
- Applications: Abdominal, Small Parts, Basic cardiol