

BACTERIAL FILTER - MOUTHPIECE for Customed, Jaeger



Product Code	33886
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	55890
CND	R049001
EAN/UPC	8023279338867

Description

BACTERIAL FILTER - MOUTHPIECE

Size - cm Ø int x ext: 2.95x3.3
Other brands: Customed, Jaeger

Disposable anatomical filter with integrated mouthpiece to prevent bacterial and viral contamination.
Easy to use, direct connection to spirometer on one side, to patient on the other side.
Suitable for most of spirometers and devices. Class II A. Individually packed in transparent bag.

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR.

Made in Italy.